

Urine Collection for Chlamydia & Neisseria by PCR 衣原體和淋病聚合酶鏈鎖反應測試

Preparation and Instruction for Investigation:

1. Collect first catch urine into a clean polypropylene container without preservatives to 1/3 full (about 30mL).
2. Seal the container and label appropriately.
3. Please send specimen together with request form to the Clinical Laboratory. If specimen cannot be delivered immediately, please keep in a cool place.

測試準備及須知:

1. 將前段尿液留於已消毒的聚丙烯容器至 1/3 瓶滿 (約 30mL)。
2. 在該容器適當位置貼上病人標籤。
3. 請把樣本連同檢驗單送往化驗室。若無法馬上送檢時，請置於陰涼處。